
norditalia

Nebulizing System Neboplus

Focused Therapy - Fast - Therapeutically Effective

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NEBOPLUS Ampul with particle selector

Particle dimension has a direct effect on drug deposit in respiratory airways.

Norms EN 13544-1 states that:

- particles with a diameter between 0,5 and 3 micron reach Lower respiratory airways fill pulmonary alveolus;
- particles with a diameter between 2 and 6 micron reach Medium respiratory airways fill tracheobronchial tree ;
- particles with a diameter over 5 micron will settle on High respiratory airways (Nose, Larynx and Trachea).

According to the above mentioned statement we created a nebulizing system adjustable according to therapy that patient needs to follow.

Neboplus Ampul is provided of multiple particle selectors that can provide a differentiation in the medicine nebulized:

- 1) for therapy aimed to act on High respiratory airways use the ampul without particle selector ;
- 2) for therapy aimed to act on all respiratory airways use the blue particle selector;
- 3) for therapy aimed to act on Lower /deep respiratory airways use the white particle selector;

Thanks to this system it is possible to increase the percentage of medicine directed towards those areas which needs to be treated, improving the effectiveness and quickness of the treatment.

Major diseases of superior tree -> Sinusitis, pharyngeal-tonsillitis, laryngitis, allergies treatments.

Major diseases of medium tree -> Bronchitis, tracheobronchial diseases, allergies treatments.

Major diseases of lower tree -> Bronchial asthma, bronchiolitis, TBC, pneumonia, alveolar disease, treatment of allergies.

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
High Therapy No selector

1 High Therapy
No selector

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
All round Therapy Blue particle selector

2 All round Therapy
Blue particle selector

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

Deep therapy White particle selector

3 Deep therapy
White particle selector

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The therapeutic result of aerosol therapy is strictly depending on the precision and to manufacturing quality of the nebulizing ampul.

The high manufacturing precision combined with the accurately studied flows dynamics, enables to reach excellent performances both in terms of nebulization speed and of therapeutic effectiveness.

The central air duct provides an additional flow that increases the nebulization speed and support the respiratory flow of patient.

The particle Selector filters the particles blocking those which have bigger size. The manufacturing precision of pixper together with the particles selector are warranty of optimal MMAD.

Each batch of production is manufactured with first choice materials and is controlled with a laser diffraction system in order to verify performances compliance.

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
The Neboplus ampul is provided of a flow regulator which enables the patient to set the nebulization speed :


- Open flow regulator -> standard speed
- Closed flow regulator -> minimum speed.

The lower speed together with an appropriate respiratory rhythm enable a higher deposition on respiratory airways, with a reduction of medicine waste.

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Neboplus ampul can work also in a reclined position, improving the therapy for those patients who are obliged to stay in bed.

It is manufactured with medical materials, biocompatible and non-toxic.

Capacity 8 ml. Residual: less than 0,5 ml

Compatible with oilish drug
Graduated scale

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Mouthpiece
Boccaglio

Tube
Tubo

All the accessory available and connectable to Neboplus ampul are manufactured in non-toxic and biocompatible materials, they have an ergonomic design for a comfortable use.

Masks are available for adults and for children. They are manufactured with a soft material and their edges are accurately molded to avoid skin irritation.

Neboplus ampul and accessory are easily washable and sterilizable.

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
LASER DIFFRACTION SYSTEM - SISTEMA A DIFFRAZIONE LASER

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Norditalia Group S.r.l.