

HOSPITAL | HOME | OUTSIDE



(A) GPRS transmission




Cell tower

DATA CENTRE



DISPATCH CENTER | EMERGENCY ROOM | CATLAB | CARDIOLOGIST



HOSPITAL | HOME



(B) Bluetooth transmission



HOME



(C) USB transmission



REALTIME MONITORING SMARTPHONE

Contact patient or hospital



(E) COMPUTER SOFTWARE FOR DATA COLLECTION, ANALYSIS AND REPORTING

(D) Software installed on a PC




(F) Dyna-Vision SwitchBoard:
Start Real-Time (bedside) Monitoring by selecting a DV unit from the switchboard

Dyna-Vision Serial Number	Status	Event Date	Event Time	Event Type	Event
210819180122	Online	01/10/2010	14:50	Alarm	Substrate
210819180122	Online	01/10/2010	14:46	Alarm	Substrate
210819180125	Online	01/10/2010	14:37	Notification	Unit is Offline
210819180122	Online	01/10/2010	14:38	Alarm	Substrate
210819180125	Online	01/10/2010	14:32	Notification	Unit is Offline
210819180122	Online	01/10/2010	14:28	Alarm	Substrate
210819180122	Online	01/10/2010	14:13	Alarm	Substrate
210819180122	Online	01/10/2010	14:14	Alarm	Substrate
210819180122	Online	01/10/2010	14:10	Notification	Unit is Offline
210819180125	Online	01/10/2010	14:05	Notification	Unit is Offline
210819180122	Online	01/10/2010	14:06	Alarm	Substrate
210819180125	Online	01/10/2010	14:05	Notification	Unit is Offline

(G) Dyna-Vision Event Dashboard:
Receive and view cardiopulmonary events directly on the desktop

(D) Software installed on a PC



(H) Resting ECG analysis
12-channel ECG analysis



(J) Holter Analysis
3 and 5 lead wire ECG holter analysis