

Nursing aid

**SAMARIT**  
safe patient moving & handling

# SEVERO

Electric pill grinder

From pill  
to powder  
in only a few  
seconds!



Quick – Easy – Safe

[www.samarit.com](http://www.samarit.com)

CE

# SEVERO

Electric pill grinder

- ✓ No more exposure to medication dust
- ✓ No more spills and cross contamination
- ✓ Prevention of RSI
- ✓ Fast and easy to use
- ✓ Near silent crushing
- ✓ Portable design
- ✓ Ideal for administering via feeding tubes



## From pill to powder...

**Many people have difficulties swallowing pills. Grinding pills and adding liquid facilitates the intake of medication.**

Manually grinding medication is time consuming work, which frequently leads to complaints of physical pain. There is also a risk of inhaling aerosols from the medication being transformed.

### Severo – quick and easy

All of these problems are a thing of the past with the Severo Pill Grinder. It's very simple. Put the medication in the SeveroCup, then place them into the grinder, close the door and within seconds the medication is available in a powder-form. The medica-

tion can also be mixed with a liquid using the same SeveroCup.

### Safe and cost effective

Did you know that grinding a Tylenol uses a force of over 21 kg (46.2 lbs). Electrical grinding is time saving and prevents complaints of physical pain. This means more time for the patients, which is an important profit in a time of increasing patient numbers and tighter economy. Grinding in a closed compartment reduces the risk of

inhaling medication powder. Save with the Severo.

Nowadays the Severo pill grinder is used in hospitals and health care facilities worldwide.

## ...in only a few seconds!

'A must-have in every health care facility!'

**SAMARIT**  
safe patient moving & handling

**Samarit Medical AG**

Gewerbestr. 12

CH-8132 Egg, Schweiz

Tel.: +41 44 918 10 11 Fax: +41 44 918 24 39

e-mail: info@samarit.com - www.samarit.com



CE

